



PROCESS SYSTEMS INTERNATIONAL, INC.

LIGO-E970192-00-V

FEDERAL EXPRESS

TO : LIGO PROJECT
MAIL STOP 18-34
CALTECH
PASADENA, CA 91125

DATE : 12/23/97
TRANS. NO. : CT148
PROJECT NO.: V59049

ATTN : LINDA TURNER

SENT BY : RICH BAGLEY

THE FOLLOWING DRAWINGS
 DOCUMENTS
 SPECIFICATIONS

ARE ATTACHED
 SENT SEPARATELY

Document No.	Rev.	Title	Dwg.Size	Sheets
01 02 03 V049-2-186	0	THIS SPEC WAS NOT PREVIOUSLY PART OF FDR. IT REPLACES V049-2-137 WHICH IS NOW VOID PROC RGA FIELD CALIBRATION LIGO-E970192-00-V	A	22

CC : JOHN WORDEN - WA SITE

ISSUED FOR: APPROVAL PRELIMINARY
 REVIEW CONSTRUCTION
 INFORMATION DISTRIBUTION

ENCL: PRINTS REPRODUCIBLES

PAGE 1 OF 1

Title: PROCEDURE FOR RGA FIELD CALIBRATION ON AN ISOLATABLE SECTION

PROCEDURE FOR RGA FIELD CALIBRATION FOR AN ISOLATABLE SECTION
LIGO VACUUM EQUIPMENT

Hanford, Washington and Livingston, Louisiana

JOB NO. V59049

PREPARED BY:

John Flier

QUALITY ASSURANCE:

Robert L. Thamm

TECHNICAL DIRECTOR:

D. C. McWilliams

PROJECT MANAGER:

Richard Bagby

DEC 23 1997

REV LTR	BY-DATE	APPD. DATE	DESCRIPTION OF CHANGE
0	JBF-12/19	REB 12/22/97	INITIAL RELEASE DEO 0581
PROCESS SYSTEMS INTERNATIONAL, INC.			SPECIFICATION
INITIAL APPROVALS	PREPARED DATE	Approved DATE	Number: V049-2-186 A LIGO-E970192-00-V
	J.F. 12/19/97	REB 12/22/97	Rev.0

OUTLINE

1.0 PURPOSE

2.0 GENERAL

3.0 REFERENCE DOCUMENTS

4.0 RESPONSIBILITY

5.0 PARTIAL PRESSURE MEASUREMENT IN AN ISOLATABLE SECTION

5.1 Option 2: Calibrated Mixed Gas Leak with Fixed Orifice

6.0 CALIBRATION METHOD

6.1 Calibration at 1×10^{-7} Torr

7.0 SETUP AND RGA CONDITION

7.1 Setup

7.2 Pump-down of Calibration Chamber

7.3 Bakeout

8.0 CALIBRATION PROCEDURES

8.1 RGA Mass Scale Tuning, Ion Source Setup, and Detector Parameter File Setup

8.2 RGA and Stabil-Ion Gauge Bakeout and Soak

8.3 RGA and Stabil-Ion Gauge Cooldown

8.4 RGA and Stabil-Ion Gauge Degas

8.5 Background (Baseline) Scan

8.5.1 Bargraph Scan

8.5.2 Pressure Readings

8.5.3 Analog Scan

8.6 Mixed Gas Calibration

8.7 Shutting Down the Calibration System

SPECIFICATION

Number: V049-2-186

A

Rev.0

9.0 CALIBRATION OF RGA

LIST OF DATA SHEETS

- 1.0 RGA Calibration Chamber Pumpdown Log**
- 2.0 RGA Calibration Chamber Bakeout Temperature Log**
- 3.0 RGA Ion Source Settings Sheet**
- 4.0 RGA Scan Parameter File Settings**
- 5.0 RGA Calibration Chamber Partial Pressure Data Sheet**
- 6.0 RGA Computer Data File Log**

SPECIFICATION

Number: V049-2-186

A

Rev.0

1.0 PURPOSE

The purpose of this procedure is to define the steps necessary to calibrate an RGA for determining partial pressures of the various gas species in the vacuum on a large volume (isolatable section).

2.0 GENERAL

This procedure is generally applicable for any RGA, but specific reference will be made to the Balzers Quadstar software and for the Balzers QMS 200 PRISMA RGA.

Data acquisition and control of the RGA is done with a PC through the RS-232 interface using the software provided with the RGA.

Software should have been loaded on the computer to be used.

3.0 REFERENCE DOCUMENTS

Balzers QUADSTAR 421 SOFTWARE MANUAL
PSI Spec. # V049-2-113, V049-2-114, V049-2-115

4.0 RESPONSIBILITY

The procedure is applicable to PSI Personnel.

SPECIFICATION

Number: V049-2-186
A

Rev.0

5.0 Partial Pressure Measurement in an Isolatable Section

5.1 OPTION 2: Calibrated Mixed Gas Leak with Fixed Orifice

An absolute calibration is required for the RGA to measure actual partial pressures of the residual gasses.

Outline of the calibration method

A mixed gas leak shall be calibrated against a NIST traceable standard for N₂. The leak should consist of H₂, N₂, Ar, and Xe. The calibration values from the vendor will be used without any corrections.

The calibrated Stabil-Ion gauge shall be mounted on the calibration chamber with a 1½-inch all-metal valve (C= 40 l/s).

The RGA shall be calibrated on the calibration chamber with a 1½-inch all-metal valve (C= 40 l/s).

Sensitivity shall be determined for the following gasses against the calibrated mixed leak with the fixed orifice: H₂, CH₄, H₂O, N₂, Ar, CO₂ and Xe. Adjustments using published ionization and transmission efficiency factors relative to nitrogen will be made for CH₄, H₂O, and CO₂.

Calibration shall be done in the 1x10⁻⁷ Torr range due to orifice size. Background pressure is expected to be about 10⁻⁸ Torr in the calibration chamber after baking pumping through the fixed orifice. The mixed gas leak of 6x10⁻⁶ Torr-L/s with about 1% Xe will be used to check RGA sensitivity in the 10⁻⁸ Torr Range at the higher AMU (Xe). Pump speed for a ¼-inch diameter orifice is about 3.7 l/s for N₂.

A data point shall be taken with the RGA and total pressure gauge to check the sensitivity for H₂O prior to bake out when the chamber is wet.

The chamber and instruments shall be baked to eliminate moisture. The baking will occur when the test chamber is attached to the BSC's 2½-inch RGA port (C=118 l/s). This port is located off the V049-4-045 or V049-4-046 manifold which is attached to the V049-4-142P1 flange.

Once calibrated, the RGA is then rebaked along with the system. The RGA is only exposed to the system when the system has reached high vacuum.

SPECIFICATION

Number: V049-2-186

A

Rev.0

6.0 CALIBRATION METHOD

The calibration method used in the field station commissioning is Option 2 from the Doc. V049-PL-485 titled "Calibrated Mixed Leak + Orifice". This is the preferred method from LIGO due to its high level of repeatability. This method of calibration assumes a linear relationship between the sensitivities at different pressure ranges (10^{-6} to 10^{-9}), as long as the electron multiplier is turned on.

6.1 Calibration at 1×10^{-7} Torr

Determine sensitivities for the selected gas species by operating the 250 l/s turbomolecular pump with the fixed orifice and the calibrated mixed gas leak open. The fixed orifice will determine the pressure profile in the test chamber. This can be compared to the measured Stabil-ion gauge reading. RGA sensitivities and partial pressures can be calculated from this known pressure profile.

7.0 SETUP AND RGA CONDITION

See attached drawing on page 8 for reference.

7.1 Setup

RGA

The RGA will be located on the calibration vessel along with a hot ion gauge. If the RGA head is valved off to the calibration vessel with an angle valve then a 1½-inch all-metal valve should be used.

RGA should be located 90° from the hot ion gauge if the RGA is in the line of sight. The preferred orientation is horizontal if the RGA is to be turned on while warm. This will prevent hot air from convecting onto the electronics unit when it is mounted on the sensor head.

Connect the RGA electronics package to the RGA detector head and the communications cable to the computer.

Mixed Gas calibrated leak

The mixed gas calibrated leak should be supplied into the vacuum space at the far end of the calibration chamber. This will allow the leak to distribute into a uniform molecular flow and to allow the proper pressure profile to be established.

SPECIFICATION

Number: V049-2-186

A

Rev.0

Title: PROCEDURE FOR RGA FIELD CALIBRATION ON AN ISOLATABLE SECTION

Hot Ion Gauge

The hot ion gauge used will be a Model 360 Granville-Phillips Stabil-Ion Gauge.

Data acquisition

Data will be taken by two methods: A complete 1-200 AMU Scan in BARGRAPH mode and one in ANALOG mode. Both these scans shall be made with the Electron Multiplier on.

7.2 Pump-down of Calibration Chamber

Rough the chamber down with the Auxiliary Turbo Cart with the turbomolecular pump operating. Start the Main 250 l/s turbomolecular pump.

Record the pressure vs. time on the data sheet. This is done to have a history on roughing the calibration chamber. The operator can recognize a problem for subsequent calibrations in the event a valve was left partially open or a large leak develops.

Verify that the test chamber is leak tight ($<1 \times 10^{-9}$ Torr-l/s) prior to bakeout.


7.3 Bakeout

Prior to calibration, the RGA, Stabil-Ion Gauge, and calibration chamber must be baked. For the Balzers PRISMA RGA, the detector head can be baked to 200°C with the electronics package removed. Requirements for bakeout are that warm-up of the RGA head shall not occur until the pressure is below 10^{-4} Torr to prevent the bake-on of contaminants. Bake for 24-48 hours. Start calibration procedure when all components are at room temperature.

SPECIFICATION	
Number: V049-2-186 A	Rev.0

**PSI Proposed
RGA Calibration
Device Utilizing
Option 2 from
Doc. V049-PL-485**

Calibration at 10^{-7}
pressure range with EM on



SPECIFICATION	
Number: V049-2-186 A	Rev.0

8.0 CALIBRATION PROCEDURES

Summary of Procedures

- 8.1 RGA Mass Scale Tuning, Ion Source Setup, and Detector Parameter Files Setup
- 8.2 RGA and Stabil-Ion Gauge Bakeout and Soak
- 8.3 RGA and Stabil-Ion Gauge Cooldown
- 8.4 RGA and Stabil-Ion Gauge Degas
- 8.5 Background (Baseline) Scan
 - 8.5.1 Bargraph Scan
 - 8.5.2 Pressure Readings
 - 8.5.3 Analog Scan
- 8.6 Mixed Gas Calibration
- 8.7 Shutting Down the Calibration System

8.1 RGA Mass Scale Tuning, Ion Source Setup, and Detector Parameter Files Setup

MASS SCALE TUNING

Mass scale tuning should be done only if the RGA has not been used for a long time, if a filament has been replaced, or the second filament chosen for use (there are two filaments available on each head). Mass scale tuning allows one to calibrate the detectors mass scaling against a known source.

Open the calibrated mixed gas leak.

Select Program Icon "TUNE-UP", submenu "Tune", and selection "QMS200 Tune mass scale".

To align the RGA mass scale with the peaks from the gas source, two parameters need to be adjusted: the offset and slope.

The offset shifts the mass scale axis left or right.

The slope adjust (stretches the axis) spacing between AMU tick marks relative to the actual calibration peaks. It may not always be possible to align the mass scale axis to the actual peaks perfectly. The calibration gas gives peaks for Hydrogen, Nitrogen, Argon, and Xenon.

Open the calibrated leak and align the Mass scale axis with the peaks, going back and forth between the higher AMU and lower AMU peaks.

Table 8.1 defines atomic mass number for each gas species of interest.

SPECIFICATION	
Number: V049-2-186 A	Rev.0

Title: PROCEDURE FOR RGA FIELD CALIBRATION ON AN ISOLATABLE SECTION

Table 8.1

Species	Mass No.
H ₂	2.0
CH ₄	16.0
N ₂	28.0
Ar	40.0
CO ₂	44.0
Xe	129
	131
	132
	134
	136

ION SOURCE

Ion source setup is located under Program Icon "TUNE-UP", submenu "Tune", and selection "Ion Source". The Ion source settings allow you to set the filament current, filament (#1 or #2), and the filament protection current. It also allows for setting of voltages, which are factory set and should not be changed without consulting the proper personnel.

Type	CH-TRON	IS-TYPE:	HS-THOR.
Channel	0 ENABLE		
Detector		Amplifier	RF-Polarity
Type	CH-TRON	Range	inverse
SEM Volt.	<< 1700>>	Offset	ON
Mass		Ion Source	IS-Voltages
Mode	SCAN-N	Filament #	[V]
First	0.00	IS-Set	IonRef
Width	6		138
Speed	5s	IS-Emission	Cathode
Resolution	25	Emiss [mA]	90.0
Threshold	-	Protect [A]	Focus
			9.38
			Field Axis
			5.75
			Extract
			12
			Fil.Prot.
			Thresh.
			[mbar]
			ON below
			OFF above

SPECIFICATION

Number: V049-2-186

A

Rev.0

Title: PROCEDURE FOR RGA FIELD CALIBRATION ON AN ISOLATABLE SECTION

PARAMETER FILE SETUP

ANALOG SCAN PARAMETER FILE:
LIGO200.SAP

BARGRAPH SCAN PARAMETER FILE:
LIGO200.SBP

Load Ch 00	CH-0
State	ENABLE
Det. Type	CH-TRON
Mass Mode	SCAN-F
First Mass	0.00

Load Ch 00	CH-0
State	ENABLE
Det. Type	CH-TRON
Mass Mode	SCAN-F
First Mass	0.00

Detector

SEM Voltage	1700
-------------	------

Detector

SEM Voltage	1700
-------------	------

Mass

Speed	5 s
Width	200
Resolution	25

Mass

Speed	5 s
Width	200
Resolution	25
Threshold	1E-15

Amplifier

Amp. Mode	AUTO
Amp. Range	----
Range-L	----
Pause - Cal	1.0
Offset	ON

Amplifier

Amp. Mode	AUTO
Amp. Range	----
Range-L	----
Pause - Cal	1.0
Offset	ON

OUTPUT: User discretion
DISPLAY: User discretion

SPECIFICATION	
Number: V049-2-186 A	Rev.0

8.2 RGA and Stabil-Ion Gauge Bakeout and Soak

RGA Bakeout

The RGA head needs to be baked along with the isolatable section. The RGA shall be baked at a temperature of 160°C or higher.

Stabil-Ion Gauge Bakeout

The Stabil-ion gauge will be on during the bake-out to monitor the pressure. The Stabil-ion gauge will be baked at a temperature of 160°C or higher (same setpoint as RGA).

8.3 RGA and Stabil-Ion Gauge Cooldown

These components should soak for 2 hours longer than the chamber soak time. The ramp down should be staggered in order to maintain the gauge temperature 25°C above the chamber temperature.

8.4 RGA and Stabil-Ion Gauge Degas

Stabil-Ion Gauge Degas

Degas the hot ion gauge for 4 min at the end of the bake soak cycle, before the start of cool down.

RGA Degas

When the vacuum vessel temperature is below 70°C, the heat on the RGA head can be turned off and the RGA head allowed to cool. When the RGA head temperature drops below 100°C, the electronics unit can be mounted onto the head and the RGA can be turned on and degassed.

Execute the program icon **TUNE-UP** and select the “**connect**” option under “**Comm**” menu to connect to the RGA. Once communications is established, select “**degas**” option under the “**setup**” menu.

Select the following degas settings:

Degas Control QMA200

Degas	
Filament #	Fil 1 or Fil 2
Current [mA]	10.0
Protect [A]	3.50
Time [min]	4

SPECIFICATION

Number: V049-2-186
A

Rev.0

8.5 Background (Baseline) Scan

After cooldown of the chamber and instrumentation, the background pressure for the chamber can be taken.

8.5.1 Bargraph Scan

Start scanning in the bargraph mode after degas of the RGA when the chamber is still warm.

Execute QUADSTAR Software "Measurement" icon and turn on filament and electron multiplier if not already on.

Select "Scan" and submenu "Bargraph".

Load the parameter file: LIGO200.sbp.

Verify settings according to Bargraph Scan settings sheet in Section 9. Enter any changes on the data sheet.

Save data to **WLERGA_1.sbc** where **WLE** stands for **W** (ashington), **L** (eft), and **E** (nd). Enter RGA serial #, SEM voltage, etc. in File Info section in the Save menu.

Verify that the Stabil-Ion gauge is on.

Continue to take Bargraph scans with the Stabil-Ion gauge on until chamber is at room temperature.

Wait at least 6 hours to ensure that the chamber is at equilibrium before continuing on to Sect. 8.5.2.

8.5.2 Pressure Readings

Take a reading of the Stabil-ion gauge at filament current of 4.0 mA setting

Wait 10 minutes and take a reading at a filament current of 0.1 mA setting.

Turn on the RGA filament and let the system stabilize for 10 minutes.

Record the total pressure at the 0.1 mA setting with the RGA operating.

The pressure in the calibration chamber should be hydrogen dominate. This means that a correction is needed for the gauge, since it is calibrated against nitrogen. Divide the Stabil-ion gauge pressure reading by the ionization efficiency for hydrogen. If the pressure in the chamber is dominated by another gas then multiply that gas's ionization efficiency by the pressure reading. See "RGA Calibration Chamber Partial Pressure Data Sheet" on Page 22 for ionization efficiency values.

SPECIFICATION

Number: V049-2-186

A

Rev.0

Title: PROCEDURE FOR RGA FIELD CALIBRATION ON AN ISOLATABLE SECTION

8.5.3 Analog Scan

Execute QUADSTAR Software "Measurement" icon and turn on filament and electron multiplier if not already on.

Select "Scan" and submenu "Analog".

Load the parameter file: LIGO200.sap

Verify settings according to Analog Scan settings sheet in Section 9, Page 21. Enter any changes on the data sheet.

Take Analog scans with the Stabil-Ion gauge on (2 minimum).

Save data to **WLERGA_1.sac** where **WLE** stands for **W**(ashington), **L**(eft), and **E**(nd). Enter RGA serial #, SEM voltage, etc. in File Info section in the Save menu. The file name should be the same as the name used in the bargraph mode.

Turn off the Stabil-Ion gauge.

Continue to take Analog scans with the Stabil-Ion gauge off (2 minimum).

The baseline pressure is expected to be predominantly H₂. With an approximate surface area of 3000 cm², the gasload should be approximately 3x10⁻⁸ Torr-L/s. With a pumping speed of 13.8 l/s (H₂ corrected) the pressure should be mid 10⁻⁹ to low 10⁻⁸ Torr range.

8.6 Mixed Gas Calibration

Open the mixed gas calibrated leak after the baseline is established.

Repeat Section 8.5.3 with the calibrated leak open for the Analog scan only.

The time constant for the calibrated leak for this small volume is less than 1 minute. The time to complete one scan at a 5-second dwell time is 16 minutes, therefore, by taking two scans, equilibrium is ensured.

8.7 Shutting Down the Calibration System

The RGA can now be shut down. The filaments and electron multiplier must be turned off from the Setup menu; they do not shut off automatically. All other components (turbo pumps, gauges, etc.) should be left running until all calculations are complete and deemed valid. Once it is complete the 2½-inch valve above the 250 l/s turbo can be closed and the 2½-inch valve to the BSC can be opened, if and only if the BSC is under a high vacuum.

SPECIFICATION

Number: V049-2-186
A

Rev.0

9.0 CALIBRATION OF RGA

Sensitivities at other AMU

To determine the pressures for species other than H₂, N₂, Ar, and Xe, previously determined sensitivities for that particular gas specie would be used. This utilizes published ionization efficiencies (E_{amu}) and transmission efficiency factors (F_{amu}). These values will be used without any corrections.

Actual orifice speed at other AMU

Determination of the orifice speed for other AMUs is dependent on the knowledge of the pumping speed for air as a reference point. The governing equation for an orifice size of ¼-inch diameter is

$$S_{orifice_amu} = 11.6A \sqrt{\frac{28.78}{amu}} = 3.673 \sqrt{\frac{28.78}{amu}} \quad \text{where the units for } A \text{ is cm}^2.$$

Calculation of Partial Pressures

- E_{amu} : Ionization efficiency for a specific AMU (-)
- F_{amu} : Transmission efficiency for a specific AMU (-)
- I_{amu} : Ion current at specific AMU of the background (Amp)
- I_{leak} : Ion current at specific AMU with calibrated leak opened (Amp)
- S_{p_amu} : Pressure sensitivity of specific AMU (Torr/Amp)
- $S_{orifice_amu}$: Pumping speed of specific AMU of that component through the orifice (l/s)
- Q_{amu_leak} : Calibration load for this AMU with calibrated leak opened (Torr-l/s)
- P_{amu} : Partial pressure of specific AMU

Sensitivities for the various gas species from the calibrated leak can be calculated from

$$S_{p_amu} = \frac{Q_{amu_leak}}{S_{orifice_amu} (I_{leak} - I_{amu})}$$

For species that are calibrated directly using the calibrated leak, the specie's system pressure is computed by

SPECIFICATION	
Number: V049-2-186 A	Rev.0

Title: PROCEDURE FOR RGA FIELD CALIBRATION ON AN ISOLATABLE SECTION

$$P_{amu} = I_{amu} \cdot S_{P_{amu}}$$

Species that do not have a direct calibration will be calculated using the values for nitrogen with the ionization and transmission efficiency correction factors. The equation for the partial pressure of other AMUs is

$$P_{amu} = \frac{I_{amu} \cdot F_{amu} \cdot S_{PN_2}}{E_{amu}} \cdot \sqrt{\frac{28}{amu}}$$

For the "OTHER" category, we will add all relevant peaks and use the values for nitrogen to calculate partial pressure.

Refer to PSI Specifications V049-2-113, V049-2-114, and V049-2-115 for cleanliness level criteria.

SPECIFICATION

Number: V049-2-186
A

Rev.0

Title: PROCEDURE FOR RGA FIELD CALIBRATION ON AN ISOLATABLE SECTION

TITLE	RGA ION SOURCE SETTINGS SHEET
DATE:	
TIME:	
TEST I.D.: e.g. WBSC1_1	
PSI TEST ENGINEER:	

Type	CH-TRON	IS-TYPE:	HS-THOR.
Channel	0 ENABLE		
Detector		Amplifier	
Type	CH-TRON	Range	
SEM Volt.	<< >>	Offset	ON
Mass		Ion Source	
Mode	SCAN-N	Filament #	
First		IS-Set	
Width			
Speed		IS-Emission	
Resolution		Emiss [mA]	
Threshold	-	Protect [A]	3.5
		RF-Polarity	inverse
		IS-Voltages	[V]
		IonRef	138
		Cathode	90.0
		Focus	9.38
		Field Axis	5.75
		Extract	12
		Fil.Prot.	Thresh.
			[mbar]
		ON below	
		OFF above	

SPECIFICATION

Number: V049-2-186

A

Rev.0

Title: PROCEDURE FOR RGA FIELD CALIBRATION ON AN ISOLATABLE SECTION

TITLE	RGA SCAN PARAMETER FILE SETTINGS
DATE:	
TIME:	
TEST I.D.: e.g. WBSC1_1	
PSI TEST ENGINEER:	

PARAMETER FILE: LIGO200.SBP

PARAMETER FILE: LIGO200.SAP

Actual

Actual

Load-Ch:00	CH-F	
State	ENABLE	
Det. Type	CH-TRON	
Mass Mode	SCAN-F	
First Mass	0.00	

Load-Ch:00	CH-F	
State	ENABLE	
Det. Type	CH-TRON	
Mass Mode	SCAN-F	
First Mass	0.00	

Detector

Detector

SEM Voltage	1700	
-------------	------	--

SEM Voltage	1700	
-------------	------	--

Mass

Mass

Speed	5 s	
Width	200	
Resolution	25	
Threshold	1E-15	

Speed	5 s	
Width	200	
Resolution	25	

Amplifier

Amplifier

Amp. Mode	AUTO	
Amp. Range	----	
Range-L	----	
Pause Cal	1.0	
Offset	ON	

Amp. Mode	AUTO	
Amp. Range	----	
Range-L	----	
Pause Cal	1.0	
Offset	ON	

OUTPUT: User discretion
DISPLAY: User discretion

SPECIFICATION

Number: V049-2-186
A

Rev.0

Title: PROCEDURE FOR RGA FIELD CALIBRATION ON AN ISOLATABLE SECTION

TITLE	RGA CAL. CHAMBER PARTIAL PRESSURE DATA SHEET
DATE:	
TIME:	
TEST I.D.: e.g. WBSC1_1	
PSI TEST ENGINEER:	

AMU	I_{amu} (Amp)	Q_{amu_leak} Leak rate (Torr-L/s)	F_{amu} Transmission Factor wrt N_2	E_{amu} Ionization efficiency wrt N_2^*	I_{leak} (Amp)	S_{p_amu} (Torr/A)	P_{amu} (Torr)
H ₂		4.8x10 ⁻⁶		0.46			
He				0.18			
12			0.42				
14			0.5				
15			0.54				
CH ₄			0.57	1.60			
17			0.6				
H ₂ O			0.64	1.12			
19			0.67				
26			0.71				
28		9.5x10 ⁻⁷		1.00			
32			1.14	1.01			
38			1.36				
40		9.4x10 ⁻⁸		1.29			
43			1.53				
44			1.57	1.42			
129		2.2x10 ⁻⁸		2.87			
131		1.8x10 ⁻⁸		2.87			
132		2.2x10 ⁻⁸		2.87			
134		9.0x10 ⁻⁹		2.87			
136		8.0x10 ⁻⁹		2.87			
Other				1.00			

* Values used from Granville-Phillips for their B-A ion gauges (CH₄ from Leybold)

SPECIFICATION

Number: V049-2-186
A

Rev.0

Title: PROCEDURE FOR RGA FIELD CALIBRATION ON AN ISOLATABLE SECTION

TITLE	RGA COMPUTER DATA FILE LOG
DATE:	
TIME:	
TEST I.D.: e.g. WBSC1_1	
PSI TEST ENGINEER:	

BARGRAPH DATA FILE NAME:

ANALOG SCAN DATA FILE NAME:

OTHER DATA FILES:

PRINTOUTS OF:

1. Baseline Analog scan with gauge on.
2. Baseline Analog scan with gauge off.
3. Calibration Analog scan with gauge on.
4. Calibration Analog scan with gauge off.
5. Baseline Bargraph scan with gauge on.
6. Baseline Bargraph scan with gauge off.
7. Calibration Bargraph scan with gauge on.
8. Calibration Bargraph scan with gauge off.

SPECIFICATION

Number: V049-2-186

A

Rev.0