

Revisions to the "work in progress" v2 procedure were requested by Garilynn and Rick on 28 Nov 2016. These steps of the procedure were removed, due to the preference for a stacked approach (rather than a sub-assembled approach). They are included here as an appendix, in case

- BOLT THE SPRING PLATE TO THE MIRROR MOUNT, WITH CONTACT SURFACES BENT IN TOWARD THE MIRROR MOUNT.
- PLACE THE SPRING PLATE WITH FLAT SIDE DOWN ON A TABLE AND CONTACT SURFACES POINTING UP
- PLACE THE MIRROR ONTO THE CONTACT SURFACES OF THE SPRING PLATE, VISUALLY CENTERED WITH RESPECT TO THE INNER DIAMETER OF THE MIRROR MOUNT. 3.
- DRIVE THE SET SCREWS AND PLUNGER FORWARD INTO THE BARREL OF THE OPTIC. IT MAY BE HELPFUL TO COUNT REVOLUTIONS TO ASSIST IN THE CENTERING OF THE OPTIC, WHERE EQUAL REVOLUTIONS WILL ENACT EXACT CENTERING. ALTERNATIVELY, SET SCREW AND PLUNGER SIZE AND DEPTH MAY BE MEASURED AND USED TO ESTABLISH A TRUE CENTER.
  - THE FINAL POSITION OF THE PLUNGER THAT GIVES GOOD HOLD IS APPROXIMATELY 1/4 TURN FROM THE FIRST CONTACT. (THIS IS AN INITIAL FINDING, AND THERE MAY BE A SUPERIOR POSITION FOR THE PLUNGER)
- (IN THE CASE OF THE D1600244 PZT MOUNT ONLY) PLACE THE PZT INTO CONTACT WITH THE {MIRROR OR CONTACT ELEMENT} AND CENTER WITH RESPECT TO THE OPTIC, USING THE TECHNIQUE DESCRIBED IN THE PRECEDEING STEP.
- INSTALL OPTIC MOUNT
  - PLACE THE HARDENED SSTL BALLS (ITEM 4) IN THE BORES OF THE BODY (ITEM 1). IT MAY BE NECESSARY TO PROP UP THE BODY AT A SLIGHT INCLINE TO ENSURE THAT THE BALLS REMAIN IN THE BORES.
  - CAREFULLY BRING THE ASSEMBLY [2, 3, OR 5] WITH INTO POSITION NEXT TO THE SPACER, WITH THE MIRROR SUPPORTED UNDERNEATH BY THE SPRING PLATE. 2.
  - TILT THE ASSEMBLY [2, 3, OR 5] INTO THE SPACER, FROM FLAT TO VERTICAL, WITH THE SET SCREWS SUPPORTING THE WEIGHT OF THE MIRROR THROUGHOUT. PROCEED UNTIL THE FACE OF THE MIRROR (OR PZT, WHERE APPLICABLE) CONTACTS THE HARDENED SSTL BALLS [4]. 3.
  - WHEN THE MIRROR HAS BEEN PLACED INTO CONTACT WITH ALL THREE HARDENED SSTL BALLS [4], USE ONE HAND TO HOLD THE ASSEMBLY [2, 3, OR 5] IN PLACE. WITH THE OTHER HAND, BOLT THE MIRROR MOUNT TO THE BODY [1] WITH THE M2.5 SCREWS [ITEM 7].

LIGO CALIFORNIA INSTITUTE OF TECHNOLOGY MASSACHUSETTS INSTITUTE OF TECHNOLOGY SCALE: 1:4 PROJECTION: SHEET 2 OF 2